TREADFORTH
GLOBAL VENTURES INC.

iKEEPFIT



AirFlow Boots

User Manual

Warranty

One year limited warranty

If there are quality issues during the warranty period, we will provide free repair or replacement. Damage caused by improper use or non-standard operation is not covered under warranty. Please retain the warranty card when requesting service.

Warranty does not apply in the following situations

1. Product damage caused by improper operation.
2. Product damage caused by self-disassembly.
3. Incidental damages (such as loss of business interests, user claims, etc.).

User Warranty Card

Warranty	1 Year	Purchase Date	
User	Address		
	Name		
	Contact No.		
Reason			

Troubleshooting

If the device is not functioning as expected, please check the following. If the suggestions below do not resolve the issue, please return the product to your distributor.

Issue	Possible Cause	Solution
Unable to start	Low battery	Please charge device
	Power button not held firmly	Long press the power button
Buttons unresponsive	Liquid on surface	Wipe with a dry soft cloth
Unable to charge	USB not properly connected	Ensure USB is fully inserted
	Power source insufficient	Ensure power source outputs correctly

Maintenance

AirFlow Boots

Product Care

Please protect the touch-sensitive display screen. Clean the surface with a soft, dry cloth only. Do not use liquids or abrasive materials.

Storage

Please keep the product away from moisture, high temperatures, direct sunlight, and water. Keep out of reach of children. When not in use for an extended period, please charge it once every two months to maintain battery health.

Proper Disposal


The lithium battery used in this product is a valuable and recyclable resource.

Please remove the battery and hand it to a professional recycling organization when disposing of the product.

Dispose of expired batteries in accordance with local regulations or environmental guidelines.



Key Descriptions



Power on	Indicators stay on
Low Battery	Indicators show fast blinking
Charging	Indicators show slow blinking
Full Charge	Indicators show steady light
Lock Screen	Airbags indicators blink

Store safely, keep away from children


Getting Started

How to use

Step 1

Long press the power button to turn on.

(Default setting: M1 mode, intensity 60mmHg, 15-minute session timer.)



Step 2

Tap the timer button to adjust session time in 5-minute increments, up to 30 minutes. Tapping at 30 minutes resets to 10 minutes.

(Default: 15 minutes.)



Step 3


Tap the **"+/–"** button to adjust the massage intensity, with 5 levels: 60 mmHg → 95 mmHg → 130 mmHg → 165 mmHg → 195 mmHg. Default intensity is 60 mmHg.

Start at the lowest setting for your first session.



Step 4


Press **"</>"** to switch massage modes: M1, M2, M3. Default is M1. In M1/M3 mode, pressing **"</>"** again keeps the same mode.



Step 5


Tap the manual control to choose the massage area. Each tap cycles through:

Airbag 1 → Airbag 2 → Airbag 3 → Airbag 4 → Airbag 5 → All airbags → Airbag 1, and repeats. Hold the button to switch back to M1 mode.




Step 6

During a session, tap the power button to lock/ unlock the screen. In lock mode, the airbag and position indicators will blink.



Step 7

After your session, hold the quick deflate button to rapidly empty the airbags for easy storage.



Product Specifications

Product Name	AirFlow Boots
Model	KPF-Leg 12
Battery	4400mAh
Session Timer	15 minutes (adjustable 15-30 minutes)
Rated Voltage	7.4V
Charging	5V 2A (USB Type-C)
Rated Power	14W

Operating Environment

	Operating Environment	Storage and Transport
Temperature Humidity	5°C ~ 45°C < 80% RH 70 KPA – 106 KPA	-10°C ~ 50°C < 80% RH 70 KPA –110 KPA

Caution

| Extended sessions beyond recommended duration may reduce comfort – please use as directed

| Default auto shutdown: 15 minutes (adjustable up to 30 minutes). Restart to continue your session

| When starting, begin at the lowest intensity setting and gradually adjust to a comfortable level.

| Please use over clothing to avoid direct skin contact

| Do not stand or walk while using the device

| If you experience dizziness, headache, shortness or breath or any physical discomfort during use, stop immediately

| Please turn off device and release air from sleeves after each session, and store properly

| This device is a personal wellness massager and does not replace any other product or service

Comfort Advisory

| Do not use if the leg area has active skin irritation, recent cuts, open wounds or bruising

For best comfort, this product should not be used under the following conditions:

- Active inflammation on the calf area
- Acute skin infections
- Open wounds or recent lacerations on the calf area
- Calf area with lesions or growths
- if you prefer not to experience a compression sensation in the lower limbs

Warning

The following should not use this product:

1. People prohibited from massage by a doctor (people with blood clots, severe arterial aneursyms, acute venous aneurysms, and skin infections)
2. People with osteoporosis, spinal fractures, and muscle strains
3. People wearing pacemakers, or people with implanted electronic medical devices
4. Pregnant women and women who have just completed child birth
5. People with severe circulatory disorders around their legs
6. People with sensory impairments should not use this device without a guardian

People who have any of the following conditions should consult a doctor before using this device

1. People with malignant tumors
2. People with heart disease
3. People with loss of thermal sensation
4. People with wounds in the calf
5. People with a body temperature above 38C (with fever symptoms)
6. This device uses airbag compression. People with limb numbness or reduced sensitivity to pressure must use caution to avoid accidents, injury or discomfort.

Wellness Advisory

The following individuals should consult their wellness professional before use:

(1) Individuals advised against compression massage by a health professional

(2) Individuals with bone density concerns, spinal sensitivities, or muscle strains

(3) Individuals with implanted electronic devices such as pacemakers

(4) Pregnant individuals or those who have recently given birth

(5) Individuals with significant circulatory sensitivity in the lower limbs

(6) Individuals with sensory sensitivities should use only under the supervision of a responsible adult

The following individuals are encouraged to speak with a wellness professional before starting:

(1) Individuals with active health conditions being monitored

(2) Individuals with known cardiovascular sensitivities

(3) Individuals with reduced sensation in the lower limbs

(4) Individuals with active wound healing in the calf area

(5) Individuals with a body temperature above 38°C

(6) This device uses airbag compression. Individuals with pressure sensitivity should exercise caution and start at the lowest intensity.